



INANOIL Analyzers

Automated characterization of olive oils



Congresso SISSG 2022 "OLI E GRASSI
ALIMENTARI: INNOVAZIONE E SOSTENIBILITÀ
NELLA PRODUZIONE E NEL CONTROLLO"

Andrea Carretta

Senior Application Specialist Food and Environmental
June, 17th 2022 - Perugia



INANOIL Analyzers

COMMISSION REGULATION (EEC) No 2568/91

of 11 July 1991

on the characteristics of olive oil and olive-residue oil and on the relevant methods of analysis

(OJ L 248, 5.9.1991, p. 1)



Analizzatori INANOIL

SOLUTION 1:

Alkyl-esters and Waxes

SOLUTION 2:

Sterols and Alcohols

SOLUTION 3:

Stigmastadienes

Further automations

MOSH/MOAH (DIN EN 16995:2017-08)

Including on-line Aloxidation and Epoxidation

2/3 MCPD e GEs

AOCS Cd 29a/b/c-13

COMMISSION REGULATION (EEC) No 2568/91

of 11 July 1991

on the characteristics of olive oil and olive-residue oil and on the relevant methods of analysis

(OJ L 248, 5.9.1991, p. 1)

Change of perspective




The methods listed in regulation often involve:

- Large amount of solvents
- Need of skilled lab technicians for several hours
- Use of large quantities of consumables and glassware
- Laborious sample preparation
- Predisposition to random errors

Automated sample-prep and analysis via HPLC/GC




- Drastic reduction in volume of solvents to be used
- The need to dedicate operators limited to a few minutes
- Almost no use of consumables and glassware
- Sample preparation substantially limited to an initial dilution
- Elimination of random errors

Traditional and INANOIL workflow comparison – Alkylesters, Waxes, Stigmastadienes



Saponification

Traditional and INANOIL workflow comparison – Sterols, Alcohols



Saponification


Change of perspective

Traditional methods

HPLC/GC hyphenation

Principle



- 1 – HPLC injection (whole sample or saponified)
- 2 – Fraction(s) withdrawal
- 3 – Vial storing
- 4 – Further automated sample treatments
- 5 – Final GC injection



Optimisation results – Alkylesters and Waxes

ANNEX XX


Method for the determination of the content of waxes, fatty acid methyl esters and fatty acid ethyl esters by capillary gas chromatography



Optimisation results – Alkylesters and Waxes

ANNEX XX


Method for the determination of the content of waxes, fatty acid methyl esters and fatty acid ethyl esters by capillary gas chromatography



Reference chromatograms
comparison

Waxes

- I.S. = Lauryl arachidate
- 1 = Diterpenic esters
- 2+2' = C₄₀ esters
- 3+3' = C₄₂ esters
- 4+4' = C₄₄ esters
- 5 = C₄₆ esters
- 6 = Sterol esters and triterpene alcohols



Optimisation results – Sterols and alcohols


ANNEX V

DETERMINATION OF THE COMPOSITION AND CONTENT OF STEROLS AND TRITERPENES DIOLCOHOLS BY CAPILLARY-COLUMN GAS CHROMATOGRAPHY

Reference chromatograms comparison


Aliphatic Alcohols

Reference 2568/91



1. Eicosanol (IS)
2. Decosanol
3. Tricosanol
4. Tetracosanol
5. Pentacosanol
6. Hexacosanol
7. Heptacosanol
8. Octacosanol

INANOIL




Optimisation results – Sterols and alcohols

ANNEX V

DETERMINATION OF THE COMPOSITION AND CONTENT OF STEROLS AND TRITERPENES DIALCOHOLS BY CAPILLARY-COLUMN GAS CHROMATOGRAPHY


Reference chromatograms comparison

Sterols, triterpenyl dialcohols



1. Cholesterol
2. Cholestanol (IS)
3. Brassicasterol
- * Ergosterol
4. 24-methylencholesterol
5. Campesterol
6. Campestanol
7. Stigmasterol
8. Δ -7-campesterol
9. Δ -5,23-stigmastadienol

10. Clerosterol
11. β -sitosterol
12. Sitostanol
13. Δ -5-avenasterol
14. Δ -5,24-stigmastadienol
15. Δ -7-stigmastenol
16. Δ -7-avenasterol
17. Eritriol
18. Uvaol



Optimisation results – Stigmastadienes


ANNEX XVII

METHOD FOR THE DETERMINATION OF STIGMASTADIENES IN VEGETABLE OILS


*Reference chromatograms
comparison*

Stigmastadienes

Reference 2568/91



INANOIL



Reference Sample: Ring Test «RT62»
(Camera di commercio di ROMA)
Blend of virgin and refined olive oil

Validation approach:

- Batches of 10 replicates – Complete automation
- Evaluation of Average, CV_r%, BIAS%

Validation results – Alkylesters and Waxes

Accuracy and Precision

Validazione CERE

#	mg/kg (sommatoria)			
1	225.8	media 220.5	dev. St. 5.1	
2	216.9			
3	215.2	Valore di riferimento (RT62) Cere = 221.0		
4	215.5			
5	227.3	bias% 0.25%	CV _r % 2.32%	
6	215.1			
7	224.2			
8	224.7			
9	215.8			
10	224.0			

Validazione Alchilesteri

#	mg/kg FAEE (sommatoria)
1	36.6
2	35.6
3	36.0
4	35.5
5	35.0
6	35.2
7	35.0
8	35.0
9	35.2
10	35.0

media 35.4	dev. St. 0.53
Valore di riferimento (RT62) FAEE = 36.2	
bias% 2.18%	CV,% 1.51%


Validation results – Stigmastadienes

Accuracy and Precision

Validazione Stigma RT62

#	mg/kg
1	3.68
2	3.78
3	3.82
4	3.78
5	3.77
6	3.81
7	3.78
8	3.75
9	3.78
10	3.70

media 3.70	dev. St. 0.04
Valore di riferimento (RT62) Stigma = 3.70	
bias% 1.66%	CV% 1.16%



Validation results – Sterols and alcohols

Accuracy and Precision

Validazione Alcoli Alifatici

#	mg/kg (sommatoria)			
1	258.5	media 262.5	dev. St. 4.0	
2	265.5			
3	267.0	Valore di riferimento (RT62) Alcoli = 258.7		
4	263.2			
5	262.4	bias% 1.46%	CV _r % 1.53%	
6	259.0			
7	261.9			
8	256.1			
9	269.3			
10	262.0			

Validazione Steroli e Dialcoli Triterpenici

#	mg/kg steroli totali
1	1464
2	1383
3	1419
4	1398
5	1408
6	1428
7	1409
8	1439
9	1453
10	1420

media 1422	dev. St. 25
Valore di riferimento (RT62) Steroli = 1447	
bias% 1.73%	CV _r % 1.74%

Validation results – Sterols and alcohols

Accuracy and Precision

#	Colesterolo	Brassicasterolo	Campesterolo	Stigmasterolo	Δ -7-stigmastenolo	Eritrodiolo Uvaol	β -sitosterolo (totale)
1	0.125	0.047	4.22	1.69	0.610	5.76	92.7
2	0.139	0.045	4.21	1.76	0.537	5.40	92.7
3	0.125	0.047	4.18	1.76	0.596	5.75	92.7
4	0.138	0.044	4.15	1.77	0.566	5.42	92.8
5	0.139	0.045	4.14	1.75	0.570	5.44	92.7
6	0.128	0.056	4.20	1.75	0.585	5.86	92.7
7	0.121	0.053	4.17	1.76	0.579	5.52	92.7
8	0.131	0.057	4.24	1.86	0.553	5.96	92.6
9	0.118	0.044	3.91	1.76	0.597	5.42	93.0
10	0.125	0.047	4.20	1.71	0.573	5.75	92.7
RT62	0.136	0.045	3.95	1.81	0.552	5.368	92.5
Avg	0.129	0.048	4.16	1.76	0.577	5.628	92.7
CVr%	5.98%	10.27%	2.25%	2.47%	3.76%	3.72%	0.11%
BIAS%	5.10%	7.67%	5.31%	3.02%	4.45%	4.83%	0.28%

INANOIL improvements – Materials/additional devices

	Alkylesters and Waxes	Stigmastadienes	Sterols and Alcohols
Materials			
Generic lab glassware	NOT Needed	NOT Needed	NOT Needed
Flasks	NOT Needed	NOT Needed	NOT Needed
LC glass columns	NOT Needed	NOT Needed	NOT Needed
Test tubes	NOT Needed	NOT Needed	NOT Needed
Separating funnels	NOT Needed	NOT Needed	NOT Needed
TLC sheets			NOT Needed
TLC chambers			NOT Needed
UV Lamp			NOT Needed
Desiccator			NOT Needed
Additional devices			
Rotavapor (or equivalent)	NOT Needed	NOT Needed	NOT Needed
Nitrogen blowing device	NOT Needed	NOT Needed	NOT Needed
Muffle	NOT Needed	NOT Needed	NOT Needed

INANOIL improvements – Chemicals

Alkylesters and Waxes

Stigmastadienes

Sterols and Alcohols

Chemicals (estimated average amounts for a 10 sample batch)

Chemical	EU Method → INANOIL	% Saving - amount saved	EU Method → INANOIL	% Saving - amount saved	EU Method → INANOIL	% Saving - amount saved
<i>Activate silica</i>	150 g → not needed	100% - 150 g	150 g → not needed	100% - 150 g	-	-
<i>Organic solvents</i>	3000 ml → 250 ml	92% - 2750 ml	3500 ml → 200 ml	95% - 3200 ml	3000 ml → 200 ml	93% - 2800 ml
<i>Hydro(alcoholic) reagents for saponification</i>	-	-	4000 ml → not needed	100% - 4000 ml	same as EU Method (so far)	-

INANOIL improvements – Lab technician commitment


	Alkylesters and Waxes		Stigmastadienes		Sterols and Alcohols	
	Operator working time (estimated average for a 10 sample batch) - hours					
	EU Method → INANOIL	% Saving - amount saved	EU Method → INANOIL	% Saving - amount saved	EU Method → INANOIL	% Saving - amount saved
Overall (including GC analysis)	12 → 8	33% - 4	15→ 9.5	36% - 5.5	10 → 4	60% - 6
Sample prep only	4 → 0.5	87% - 3.5	10 → 0.5	95% - 9.5	7 → 0.5	92% - 6.5

INANOIL improvements – Final Remarks


Aimed to....

- Drastically reduce the consumption of solvents
- Limit the intervention of expert operators to a few minutes, reducing the entire sample-prep to the simple dilution of the initial sample (eventually saponified)
- Eliminate the use of glassware and consumables
- Effectively cancel the incidence of random errors

Change of




Data
reviewing



LESS
solvent exposure

... take care of analyst's health, and let him to exploit his skills
«away from the hood»

Totally automated sample-prep and analysis





Analizzatori INANOIL

PARTE 1:

Meccanismi e Cere

PARTE 2:

Stereoisomeri

PARTE 3:

Stigmastadieni

7 luglio

Aknowledgments



Fiorenzo



Carmine



Centro Analisi Biochimiche sas

srainstruments.com



Piattaforma INANOIL per la determinazione di Steroli e Alcoli

PARTE 2

SEPARAZIONE DELLE FRAZIONI DEI COMPOSTI ALCOLICI

1. OGGETTO

L'insaponificabile preparato come descritto nella parte 1 è frazionato nei diversi composti alcolici, alcoli alifatici, steroli e dialcoli triterpenici (eritrodiole e uvaolo).



2. PRINCIPIO

L'insaponificabile può essere frazionato per mezzo della cromatografia basica su strato sottile (metodo di riferimento), rivelato per poi raschiare ed estrarre le corrispondenti bande. Come metodo alternativo, la separazione può essere realizzata mediante HPLC con colonna di gel di silice e rivelatore UV raccogliendo le diverse frazioni. Gli alcoli alifatici e triterpenici nonché gli steroli e i dialcoli triterpenici sono isolati insieme.

HPLC separation

ALCOHOLIC FRACTION

STEROLIC FRACTION



Piattaforma INANOIL per la determinazione di Steroli e Alcoli



ALCOHOLIC FRACTION

STEROLIC FRACTION

ALCOHOLIC FRACTION

STEROLIC FRACTION

